

WBC manual count using hemocytometer

Prepared by
Hadeel Al Sadoun

Objectives



- To accurately count WBC in Chamber.
 - To perform reliable dilution of blood cells
 - To calculate the number of cells/ μL

Principle

- Whole blood collected in EDTA is diluted according to the type of cell count obtained.
- The diluted blood suspension is then placed in a chamber and the cell counted
- The count is multiplied by dilution factor and reported as number of cells per microlitter (μL) of whole blood

Material

- Hemocytometer with Neubauer grid.
- Cover glass
- Diluents
- Microscope.



Methodology


- Put the cover slip or glass slip on the top of grid area in the Chamber (use air tight technique)
- Dilute you sample:
 - 1: 20 for WBC count
 - 1:200 for RBC count and platelets
- Load your sample into the laoding area in the chamber
- Count the cells in the 4 large squares for WBC
- calculate the number of cells counted / μL

Sample dilution

- Dilution of whole blood sample:
 - Diluents:
 - Acetic acid (CH_3COOH)
 - Or : dis. H_2O
 - Purpose:
 - Dilute the amount of WBC , RBC to be able to count it. (NR RBC: M $4.3\text{-}6.2 \times 10^6 / \mu\text{L}$) (F: $3.8\text{-}5.5 \times 10^6 / \mu\text{L}$)
(NR WBC: $4.3\text{-}10.8 \times 10^3 / \mu\text{L}$)
 - To lyses the RBC and platelets (the diluents lyses also the WBC but takes longer time) (time factor is critical)


Methodology

- Dilution:
 - 1:20 dilution or 1:50 (ex: chronic leukemia)
 - $(1+19=20)$
 - $(50\mu\text{L of blood} + 950\mu\text{L diluent})$
- Loading the sample:



WBC count

The hemocytometer contains 2
Neubauer counting chamber →




Each chamber contains:

*4 WBC counting squares

*Each contains 16 squares

100 RBC= 10 Platelets= 1 WBC

Chose 90° lines, count only
the cells that on those lines
(ex: L-shape)
apply it to all squares for
maximum accuracy



Calculation

- Cells/ μL =
- no. of cells in 1 large square x Dilution factor

volume factor (0.1)

Dilution factor = reciprocal of dilution (20)

Volume factor = (width x length x height) = 0.1