

LAB (6):

DIALYSIS




KAU-Faculty of Science- Biochemistry department

Analytical biochemistry lab (Bioc 343) 2012

T.A Nouf Alshareef



nf.shareef@hotmail.com



Introduction:

- Dialysis is non chromatographic separation method used to separates molecules according to size through the use of semi-permeable membranes.
- dialysis process occurs at the membrane surface, depend on **conc. differences** between two solutions on both side of the membrane.
- the porous membrane selectively allows **smaller** solutes to pass while retaining larger species

It is a simple process in which:



Factors affect on dialysis:

- Temperature, viscosity, and pressure gradient across the membrane.
- Solvent, pore size, and the nature of the membrane.
(rate of dialysis is greatest in distilled water)

Advantages of lab dialysis:

- Very gentle conditions
- Easy operation
- Wide range of samples volume
- Many membrane types
- Inexpensive materials
- Disposable membranes & devices

Dialysis membrane materials:

- Synthetic or natural membranes are used for **diffusion** or **osmosis**.
- It allows the passage of small molecules but not larger ones (Filtration applications).
- **Membrane materials:** regenerated cellulose, cellulose acetate, polysulfone, polycarbonate, polyethylene, polyolefin, polypropylene, and polyvinylidene fluoride.
- **Membrane size:** wide range of pore size (**MWCO**: M.wt. cut-off). The main feature of this membrane is that it is porous.

Regenerated cellulose (Cellu-Sep®):

- is derived from cotton:
- has a symmetric pore structure allows small molecules to migrate in either direction.

Cellophane:

- is most commonly use dialysis material.
- is thin, transparent sheet made of regenerated cellulose.
- contains traces of sulphur compounds, metal ions and some enzymes.



- Dialysis is best carried out with freshly prepared tubing
(wet is very susceptible to attack by micro-organisms).
- For storing:
the tubing is best kept with a trace of benzoic acid in the solution

Dialysis membrane preparation:

- Boil the bag for 30 min in alkaline EDTA
(activation of membrane and to make the dialysis sac soft)
- After boiling, the tubing is washed with distilled water and soaked in water.


Lab practice:

Aim:

- Separation of **maltose** from mixture solution using dialysis process.
- Our sample: (solution of starch, salivary amylase and maltose)
- dialysis bag: cellophane sealed with string by knotting.
- Solvent: distilled water.


Principle:

- Starch consists of:
 - Amylose: (M.wt 50,000 g)
 - Amylopectin: (M.wt 1,000,000 g)
- neither of them will pass through dialysis membrane



- Maltose (smaller than membrane pores) —→ moves freely across out of the membrane into the solvent
- Starch (Larger molecule) —→ is retained inside the dialysis bag


The dialysis tubing is sitting in a container of distilled water.



Dialysis tubing was tied shut, so that molecules can only escape through the pores in the membrane.


Dialysis process:

- the volume of the solvent outside the bag is **greater** than inside, over time, most of maltose will leave the bag (**Osmosis**).
- The passage of maltose into or out of the bag is in the direction of decreasing concentration, therefore displaying diffusion.
- Equilibrium is achieved when the concentration of maltose **equal** inside and outside of the bag.



Chemicals & other material:

- Salivary amylase (1 ml saliva diluted to 5 ml with distilled water)
- Iodine solution (5 mmol/L in 30 g/L Potassium iodide)
- Soluble starch (20 g/L)
- Sodium chloride solution (1g/L) buffered with 0.02 mol/L Sodium phosphate pH 6.8).
- Fehling's solution.




Procedure:




1- Label 14 small test tubes as the following:

add 10 drops of **iodine solution**



(Test for starch).

1ml of **Fehling's**




(Test for maltose).

Starch + iodine → (Blue color).

Maltose + Fehling's $\xrightarrow[\text{5 min}]{\text{Heating}}$ **Red ppt**

2- Seal one end of the dialysis bag (the bottom) with a string.



3- Fill the sac with the **reaction mixture**:
(**2.5 ml** starch, **0.5 ml** buffer, **0.5 ml** amylase)

4- Quickly (at 0 min) pipette 1.5 ml of reaction mixture (0.5 ml to tube 1S) and (1 ml to tube 1M)


5- Immerse the dialysis bag in beaker contain 100 ml dis. H_2O .

6- Place the beaker on a stir plate (stir 90 min)
(Stirring help in entry of water inside the sac and
small molecules to getting out of out of the sac)


7- every 15 min, pipette 1.5 ml from the solvent
(outside) and examine each fraction for the presence
of starch or maltose.

8- At the end of experiment (90 min):
cut the dialysis sac carefully pipette 1.5 ml from solution inside, examine them for starch
and maltose.

100 ml
dis. H_2O



Test for starch: Change of the color of iodine during separation




**75 min
(inside)
yellow color
7 S**

**15 to 75 min
(outside)
yellow color
2S- 6S**

**0 min
(inside)
blue color
1S**

Test for maltose: Formation of red precipitate during the separation



75 min (inside)
red ppt
7 M

15 to 75 min (outside)
Red ppt
2 M to 6 M

0 min (inside)
red ppt
1M

Result sheet

- Make a table of time intervals (min), and observe result for each fraction with iodine and Fehling's tests.
- Note observation in table for Fehling's test as follows:
 - (+) for a few ppt,
 - (++) for more ppt,
 - (+++)for much more ppt and so on.

Tabel : 1

Time (min)	Inside tubing		Comment
	I ₂ test	Fehling's test	
0			

Tabel : 2

Time (min)	Outside tubing		Comment
	I ₂ test	Fehling's test	
15			
30			
45			
60			
75			
90			

Tabel : 3

Time (min)	Inside tubing		Comment	Outside tubing		Comment
	I ₂ test	Fehling's test		I ₂ test	Fehling's test	
90						