#5 Experiment : Wheatstone Bridge

Name	
University ID	
Section	
Lab time	Day (), time from () to ()

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

fig 2

Theory

The Wheatstone bridge is a circuit used to compare an unknown resistance with a known resistance. A Wheatstone Bridge is 4 resistors connected in the form of a square with galvanometer connected across one diagonal and a voltage source across the other, as shown in the circuit fig. When the Galvanometer indicates zero current, the bridge is balanced. under these conditions, the voltage drop from A to C must be the same as that from A to D. Similarly, the voltage drop from B to D must be the same as that from B to C.

$$I_1 R_1 = I_2 R_3 \quad \dots 1 \quad \quad I_1 R_2 = I_2 R_4 \quad \dots 2$$

By divided 1 by 2

$$\frac{R_1}{R_2} = \frac{R_3}{R_4} \quad \longrightarrow \quad \frac{R_1}{R_2} = \frac{\rho \frac{L_1}{A_x}}{\rho \frac{L_2}{A_x}}$$

$$\frac{R_1}{R_2} = \frac{L_1}{L_2}$$

Procedure

1. We connect the circuit fig 2 and make sure that the volt switch is at the lowest value and that the current switch is at the highest value.
2. turn on the battery (voltage source) and place the voltage on the specified value not exceeding 1 volt.
3. Change the keys of the resistor-box to reach the required first value in the table and then move the slider on the ruler until the galvanometer index reaches zero) the touch of the shutter that gives the reading zero in the galvanometer is divided the length of the ruler into two parts L_1 and L_2 .
4. Record of L_1 values (corresponding to the unknown resistance wire) and L_2 (corresponding to the known resistance - resistancebox-)
5. Repeat step 3 for the next resistances in the table.
6. Calculate using equation (1) R_x value at all the values of R_s in the table.
7. Find the average \bar{R}_x readings. Then calculate the area of the wire section $A_x = \pi r^2$ from the suspended where the radius of the wire section is given.
8. Calculate using the lingering (1-3) the value of the resistivity $\rho = R_x \frac{A_x}{L_x}$. then conductivity $\sigma = \frac{1}{\rho}$
9. By reference to the tables of resistivity ρ in physics books or the Internet, determine the type of wire material used in the experiment [Electrical resistivity and conductivity - Wikipedia](#) .

#5 Experiment : Wheatstone Bridge

$r = 0.0762\text{ mm} = 7.62 \times 10^{-5}\text{ m}$ $A_x = \pi r^2 = \pi \times (7.62 \times 10^{-5})^2$ $= 1.8 \times 10^{-8}\text{ m}^2$		$L_x = 20\text{ m}$	
R_s	L_1	L_2	$R_x = R_s \frac{L_1}{L_2}$
50			
60			
70			
80			
90			
$\overline{R_x} =$			
$\rho = R_x \frac{A_x}{L_x} =$		$\sigma = \frac{1}{\rho} =$	

what is the material of the wire used as unknown resistance in this experiment?

.....