




KING ABDULAZIZ UNIVERSITY

Microbial Enzymes II


Physiology of microbiology
Bio 336





Catalase is a common enzyme found in nearly all living organisms exposed to oxygen (such as bacteria, plants, and animals) which catalyzes the decomposition of hydrogen peroxide to water and oxygen.

It is a very important enzyme
in protecting the cell from oxidative damage by reactive oxygen species (ROS).



Microbial Enzymes II



1- Catalase Activity

Hydrogen peroxide (H₂O₂) is a poisonous byproduct of metabolism that can damage cells if it is not removed. All respiring organisms therefore must have some mechanism for detoxification. Catalase is one of the common methods. When hydrogen peroxide is added to a colony of catalase-producing bacteria, it is broken down into water (H₂O) and oxygen gas (O₂). that is produced can be seen as bubbles.



Microbial Enzymes II




Experiment # 1

Catalase activity test

Purpose: Identify the presence of the enzyme Catalase

- ☐ Material and Experimental Protocols:
- ☐ Observation and results:
- ☐ Conclusion:

Microbial Enzymes II



Catalase production test

Microbial Enzymes II



2- Oxidase activity

Oxidase is an oxidizing enzyme within the series of respiratory enzymes responsible for oxidative phosphorylation reactions, which can oxidize certain aromatic amines such as Phenylenediamine to form end-color products.

If bacteria is positive the paper will turn purple.

A blue- purple color produced is indicative of a positive test.

Phenylenediamine

Oxidase



perpul color



This test is used to identify microorganisms containing the enzyme cytochrome oxidase (important in the electron transport chain). It is commonly used to distinguish between oxidase negative *Enterobacteriaceae* and oxidase positive *Pseudomadaceae*. Cytochrome oxidase transfers electrons from the electron transport chain to oxygen (the final electron acceptor) and reduces it to water. In the oxidase test, artificial electron donors and acceptors are provided. When the electron donor is oxidized by cytochrome oxidase it turns a dark purple. This is considered a positive result. In the picture below the organism on the right (*Pseudomonas aeruginosa*) is oxidase positive.



Microbial Enzymes II



Experiment # 2

Oxidase activity test


Purpose: Identify the presence of the enzyme Oxidase

❑ Material and Experimental Protocols:

❑ Observation and results:

❑ Conclusion:

Microbial Enzymes II



Oxidase activity test



3- hydrolysis of Cellulose

Cellulose is a complex carbohydrate that is composed of many sugars. Some bacteria can enzymatically dissolve cellulose under aerobic conditions to glucose units. This activity is of great importance in the decomposition of plant residues in soil where microorganisms do not accumulate in the soil.

There are two types of cellulose decomposition:

- 1 - Aerobic decomposition (Debous media).
- 2 -Anaerobic degradation(Omlyanisky media)+oil.



Experiment # 3

Aerobic decomposition of cellulose

Purpose: To determine whether aerobic bacterium produces a cellulases that will hydrolyze a cellulose .

Dubos medium

- ❑ Material and Experimental Protocols:
- ❑ Observation and results:
- ❑ Conclusion:



Experiment # 4


Anaerobic decomposition of cellulose

Purpose: To determine whether anaerobic bacterium produces a cellulases that will hydrolyze a cellulose .

Omliansky's medium

- ☐ Material and Experimental Protocols:
- ☐ Observation and results:
- ☐ Conclusion:

Microbial Enzymes II



Cellulolysis



Catalase & oxidase activity

http://www.uwyo.edu/molb2021/virtual-edge/lab10/lab10_rapid.html

https://www.youtube.com/watch?v=7bkXI_iT6lQ&t=59s

<https://www.youtube.com/watch?v=PGBf7Xehcns>

<https://www.youtube.com/watch?v=brBgR589KFo>

<https://www.youtube.com/watch?v=z1HsdAQO7m4>

<https://www.youtube.com/watch?v=YSPbMdR2K70>

hydrolysis of Cellulose

<https://www.khanacademy.org/science/ap-biology/chemistry-of-life/introduction-to-biological-macromolecules/v/hydrolysis>

